Multi Drug

RAPU08A444

en

**One step test device for the simultaneous, qualitative detection of
Amphetamines, Barbiturates, Benzodiazepines, Cocaine metabolites, Methamphetamines
Morphine, Methadone, Phencyclidine, Tricyclic Antidepressants and THC metabolites
in human urine
For professional in vitro diagnostic use only.
RAPU08A444**

INTENDED USE

The **DIAsource Drugscreen-Stick RAPU08A444** is a specific arrangement of different lateral flow chromatographic immunoassays for the detection of Amphetamines, Barbiturates, Benzodiazepines, Cocaine metabolites, Methamphetamines Morphine, Methadone, Phencyclidine, Tricyclic Antidepressants and THC metabolites in human urine at the following cut-off concentrations:

TEST DEVICE Side A	SUBSTANCE	CAS-No	Cut-off Limit Value [ng / mL]	TEST DEVICE Side B	SUBSTANCE	CAS-No	Cut-off Limit Value [ng / mL]
Amphetamines (AMP)	D-Amphetamine	[51-64-9]	1000	Morphine (MOR)	Morphine	[57-27-2]	300
	3,4-Methylenedioxy-amphetamine	[4764-17-4]	2000		Codeine	[76-57-3]	300
	D,L-Amphetamine	[300-62-9]	3000		6-Monoacetylmorphine	[-]	400
	Phentermine	[122-09-8]	3000		Morphine 3-β-D-glucuronide	[20290-09-9]	1000
	L-Amphetamine	[156-34-3]	50000		Levorphanol	[77-07-6]	1500
Barbiturates (BAR)	Butabarbital	[125-40-6]	75		Hydromorphone	[466-99-9]	3000
	Butethal	[77-28-1]	100		Ethylmorphine	[-]	6000
	Phenobarbital	[50-06-6]	100		Norcodeine	[467-15-2]	6000
	Alphenal	[115-43-5]	150		Thebaine	[115-37-7]	6000
	Aprobarbital	[77-02-1]	200		Oxycodone	[76-42-6]	30000
	Amobarbital	[57-43-2]	300		Hydrocodone	[125-29-1]	50000
	Secobarbital	[76-73-3]	300	Methadone (MTD)	Methadone	[76-99-3]	300
	Pentobarbital	[76-74-4]	300		Doxylamine	[469-21-6]	50000
	Cyclopentobarbital	[76-68-6]	600	Phencyclidine (PCP)	Phencyclidine	[60124-79-0]	25
	Butalbital	[77-26-9]	2500		4-Hydroxy-Phencyclidine	[-]	12500
Benzodiazepines (BZD)	Nitrazepam	[146-22-5]	100	Tricyclic Antidepressants (TCA)	Desipramine	[50-47-5]	200
	Clobazam	[22316-47-8]	100		Imipramine	[50-49-7]	400
	Temazepam	[846-50-4]	100		Nordoxepine	[-]	1000
	Lorazepam glucuronide	[-]	150		Nortriptyline	[894-71-3]	1000
	Alprazolam	[28981-97-7]	200		Amitriptyline	[50-48-6]	1500
	Diazepam	[439-14-5]	200		Promazine	[58-40-2]	1500
	Clorazepate dipotassium	[57109-90-7]	200		Doxepine	[1668-19-5]	2000
	Norchlordiazepoxide	[-]	200		Maprotiline	[10262-69-8]	2000
	Oxazepam	[604-75-1]	300		Trimipramine	[739-71-9]	3000
	Desalkylflurazepam	[2886-65-9]	300		Clomipramine	[303-49-1]	12500
	Flunitrazepam	[1622-62-4]	400	Marijuana (THC)	Promethazine	[60-87-7]	25000
	Nordiazepam	[-]	400		11-nor-Δ ⁹ -THC-9 COOH	[-]	30
	Chlordiazepoxide HCl	[438-41-5]	750		11-nor-Δ ⁹ -THC-9 COOH	[-]	50
	Clonazepam	[1622-61-3]	800		Δ ⁸ -THC	[5957-75-5]	15000
	α-Hydroxy-alprazolam	[-]	1200		Δ ⁹ -THC	[1972-08-3]	15000
	Bromazepam	[1812-30-2]	1500		Cannabinol	[521-35-7]	20000
	Chlordiazepoxide	[58-25-3]	1500				
	Delorazepam	[2894-67-9]	1500				
	Lorazepam	[864-49-1]	1500				
	Estazolam	[29975-16-4]	2500				
	Triazolam	[28911-01-5]	2500				
	Midazolam	[59467-70-8]	12500				
Cocaine (COC)	Benzoyllecgonine	[-]	300				
	Cocaine	[53-21-4]	700				
	Cocaine ethylene [-]	[-]	12500				
	Ecgonine [-]	[-]	32000				
Methamphetamine (MET)	D-Methamphetamine [-]	[-]	1000				
	3,4-Methylenedioxy-methamphetamine	[-]	2000				
	L-Methamphetamine	[537-46-2]	8000				
	p-Hydroxy-methamphetamine	[-]	30000				
	Mephentermine	[100-92-5]	50000				

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Immerse the tips of test into specimen until chromatography is clearly visible

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REAGENTS
The **DIAsource Drugscreen-Stick** contains mouse monoclonal antibody-coupled particles and drug-protein conjugates. The control line system contains a goat anti-mouse antibody.

PANEL 25 envelops containing a defined combination of strips in a protected case, desiccant. Combination : COC / THC / BZD (BZO) / MTD (MTE) / MOP (MOR. OPN) / AMP / BAR / MET (M-AMP) / PCP / TCA

<ul style="list-style-type: none"> For Professional in vitro diagnostic use only. Do not use after the expiration date. The test device should remain in the sealed pouch until use. Do not moisten nitrocellulose membrane with urine samples. Immerse tips of test strips for at least 10-15 seconds in the urine specimen. 	<ul style="list-style-type: none"> Avoid cross-contamination of urine samples by using a new specimen collection container for each urine sample All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent The used test device should be discarded according to federal state and local regulations.
--	---

The **DiAsource Drugscreen-Stick** can be stored at room temperature or refrigerated (2 – 30 °C). The test is stable through the expiration date printed on the sealed pouch.

<input type="checkbox"/> Do not freeze.	<input type="checkbox"/> Do not use beyond the expiration date.
---	---


The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2 - 8°C for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored below -20 °C. Frozen specimens should be thawed and mixed before testing.

MATERIALS PROVIDED

- I Multi test device
- I Package insert

MATERIALS REQUIRED BUT NOT PROVIDED

- I Specimen collection container
- I Timer


Dip **DIASource Drugscreen-Stick**
into urine specimen

DIRECTIONS FOR USE


- 1 Allow the urine specimen, test device, and / or controls to reach room temperature (15 – 30 °C) prior to testing.
- 2 Bring pouch to room temperature before opening it.
- 3 Remove the Multi test from the sealed pouch and use it as soon as possible.
- 4 Immerse the tips of the Multi test vertically in the urine specimen for at least 10-15 seconds.
- 5 Do not pass the maximum line (dipping line) on the test when immersing the tips of the Multi test.
- 6 Place the test on a non-absorbent flat and clean surface. Start the timer and wait for the red lines to appear.
- 7 The result should be read at 5-8 minutes. Do not interpret after more than 10 minutes.

INTERPRETATION OF RESULTS

SIDE A

SIDE B

Negative:* Two lines appear in each window. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the concentrations of the substances detectable with the corresponding test are below the cut-off concentration (substances & cut-off concentrations see table on page 1) or that they are not present.

Positive: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the concentration of at least one of the substances detectable with the corresponding test exceeds the cut-off concentration (substances & cut-off concentrations see table on page 1).

Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test device. If the problem persists, discontinue using the lot immediately and contact distributor / manufacturer.

* **Note:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

Information about the performance characteristics are available through our technical support department, at tech.support@diasource.be

 Consult instructions for use	 Manufacturer
 Storage temperature	 Contains sufficient for n tests
 Use by	 In vitro diagnostic medical device
 Batch code	 Panel
 Catalogue number	

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